

Gmate[®] VOICE

Blood Glucose Monitoring System

Blood Glucose Meter can TALK
Audio guidance can lead you in every step




Wide Display LCD

3.5 inch large screen comforts your eyes



Speaker

The audio guidance for test results will be provided through the speaker



Data Port

Transfer test results to the PC by connecting the cable



3.5 inch Wide display



Audio guidance



Testing Alarm



Minimal Sample
Volume (0.5 μ l)

Gmate[®] VOICE

Blood Glucose Monitoring System

How Gmate[®] VOICE Works



Insert Gmate[®] Test Strip



Collect blood sample



Apply the sample



Test result will be shown in 5 seconds (with audio guidance)

Specification

Model	PG-301 (Gmate [®] VOICE) Blood Glucose Monitoring System
Test Method	Electrochemical sensor
Calibration Method	Plasma Equivalent
Sample	Capillary Whole Blood
Sample Volume	0.5 μ l
Test Time	5 seconds
Memory	500 Blood glucose test results with time and date
Result Range	20 - 600 mg/dL (1.1 - 33.3mmol/L)
Hematocrit Range	20 - 60%
Operating Temperature	10 - 40°C (50 - 104 °F)
Operating Humidity	Less than 90%
Altitude	Up to 3,048 meters (10,000 feet)
Power Source	DC 3V, Alkaline batteries (AAA), 2ea, 1.5W
Size	62.0 × 114.6 × 13.0 mm (2.44 × 4.52 × 0.51 inch)
Weight	80g (2.82 oz)



DR NTR BIOSCIENCES PRIVATE LIMITED

M63, Medical Device Park, Sultanpur, Patancheru, Telangana-502319
info@drntrbiosciences.com, sales@drntrbiosciences.com
+91 96770 90702, +91 99400 60607 | www.drntrbiosciences.com

AUTHORIZED DISTRIBUTOR